

George Pappas
CHM. 121.004
Inst. Mark Novak
Lab Report – Experiment 5
10/2/17

100
well done

Experiment Objective

The objective of this experiment was to perform multiple chemical changes to copper wire and convert it back into pure copper by the end of the experiment. At the end of the experiment, we would also calculate our percent recovery; the goal was to get 100% of the copper back.

Conclusion

We were told that this experiment was going to take up most, if not all, of the class time we had available, which it did. The part that took the most amount of time but was also most likely the reason why we had such a decent percent recovery, was vacuum filtrating the Copper (II) Oxide a second time. This allowed us to recover most, if not all, of the solid that originally was floating around in the filtrate; in fact, the filtrate was perfectly clear. However, there still may have been some solids floating around in the filtrate, this is one point of the experiment where we may have lost some of the copper. At the end of the experiment, we had .415 grams of recovered copper. Because we started off with .437 grams, our final percent recovery was 94.97% of copper.

DATA AND CALCULATIONS TABLE for a Copper Ring or Wire

| | |
|---|--|
| Mass of Beaker and Ring/Wire | 152.584 g |
| Mass of Beaker | 152.147 g |
| Mass of Ring/Wire (Theoretical Yield) | .437 g |
| Mass of Evaporating Dish (and filter paper) and Recovered Copper | 107.243 g |
| Mass of Evaporating Dish (and filter paper) | 106.828 g |
| Mass of Recovered Copper | .415 g |
| Percent Copper Recovered from Ring/Wire (Mass of Recovered Copper) x 100% (Theoretical Yield) | 94.97% 9002 |

Show calculations.

mass of Beaker and Ring/wire - mass of Beaker = mass of Ring/wire
(Theoretical Yield)

$$152.584 \text{ g} - 152.147 \text{ g} = .437 \text{ g}$$


mass of Evaporating Dish (and filter paper) and Recovered Copper - Mass of
Evaporating Dish (and filter paper) = mass of Recovered Copper

$$107.243 \text{ g} - 106.828 \text{ g} = .415 \text{ g}$$


$$\frac{(.415 \text{ g}) \times 100\%}{(.437 \text{ g})} = 94.97\%$$

QUESTIONS AND CALCULATIONS**NET IONIC EQUATIONS (Be very neat to receive full credit!)**


Copper to copper nitrate



Copper nitrate to copper hydroxide



Copper hydroxide to copper oxide



Copper oxide to copper sulfate



Copper sulfate to solid copper



Solid zinc to zinc (II) ion and hydrogen gas



QUESTIONS AND CALCULATIONS

QUESTIONS AND CALCULATIONS**CALCULATION TABLE**

Start with your theoretical yield. Use the balanced net ionic equation and stoichiometry to calculate the mass of each of the following. Show calculations.

| | |
|---|----------|
| Mass of copper in penny/ring/wire/solution (theoretical yield) $.437\text{g} \div 63.546\text{ g/mol} = .0069\text{ moles}$ | .437 g ✓ |
| Mass of copper nitrate, $\text{Cu}(\text{NO}_3)_2$ $.0069 \times 1 = .0069\text{ moles Cu}$ $.0069\text{ moles} \times 187.554\text{ g/mol} = 1.29\text{ g}$ $\text{Cu}(\text{NO}_3)_2$ | 1.29 g ✓ |
| Mass of copper hydroxide, $\text{Cu}(\text{OH})_2$ $.0069 \times 1 = .0069\text{ moles Cu}$ $.0069\text{ moles} \times 97.56\text{ g/mol} = .671\text{ g}$ $\text{Cu}(\text{OH})_2$ | .671 g ✓ |
| Mass of copper oxide, CuO $.671\text{ g} \div 97.56\text{ g/mol} = .0069\text{ moles}$ $.0069\text{ moles} \times 1 = .0069\text{ moles Cu}(\text{OH})_2$ $.0069\text{ moles} \times 79.545\text{ g/mol} = .547\text{ g CuO}$ | .547 g ✓ |
| Mass of copper sulfate, CuSO_4 $.0069\text{ moles} \times 1 = .0069\text{ moles Cu}$ $.0069\text{ moles} \times 159.602\text{ g/mol} = 1.10\text{ g}$ | 1.10 g ✓ |

QUESTIONS AND CALCULATIONS

QUESTIONS AND CALCULATIONS

List two reasons why your percent yield might be less than 100.0%

1. Some of the Copper (II) oxide went around the filter paper and into the filtrate. ✓
2. Instead of heating, the aqueous Copper (II) ion was boiling

List three reasons why your percent yield might be more than 100.0%

1. left over zinc from changing aqueous Copper (II) ion back into solid copper.
2. Not allowing the copper to completely dry. ✓
3. Drying the copper too much, this adds oxygen in the form of copper oxide

QUESTIONS AND CALCULATIONS

QUESTIONS AND CALCULATIONS

1. Digestion of a precipitate is one of the best methods to increase particle size. In this experiment, why would that be advantageous?

With increased particle size, it would be less likely to lose copper during the vacuum filtration process. It would also be easier to transfer the copper into different containers. ✓

2. A 1995 penny has a mass of 2.489 g. How much copper (in grams) is in this penny? (or What is the theoretical yield of copper from this penny?) 2.5% copper

$$2.489 \cdot .025 = 0.0622 \text{ g} \quad \checkmark$$

3. The mass of the recovered copper is 0.051 g. What is the percent yield of copper? (Note: Refer back to question #2)

$$\text{Percent Yield} = \frac{0.051 \text{ g}}{0.0622 \text{ g}} \times 100 = 81.99\% \quad \checkmark$$

